



# Electrophoresis

**DNA  
Protein  
RNA**

## UV Light Box for DNA, RNA, Protein Gel Electrophoresis




DNA Agarose Gel



### Specification

1. Wavelength: 306nm (or 254nm, 365nm)
2. Size: 50(w) x 32(D) x 15 (H) cm
3. View size: 20x20cm (or 20x35cm)
4. Light source: 15W x 6 tubes
5. Light intensity: high & low selectable
6. Fast lighting up: lighting up at once when switch on and no blink
7. Extra magnetic working area for tube rack or operating tool
8. Big size UV resistant plastic cover 30x34cm
9. Cool fan is included
10. Closed UV view : 4 sides of cover stop UV
11. Safety device: UV leakage indicator will have color change, if UV come out from plastic cover
12. Diffuse UV background: decrease UV tube image

**Special device** is more safe for user, the color of indicator will change from white to red, if detect some UV leakage from cover



### Order number:


model UV2020-1 : 254nm, view size=20x20cm  
 model UV2020-2 : 306nm, view size=20x20cm  
 model UV2020-3 : 365nm, view size=20x20cm  
 model UV2020-B : 254 /365nm, view size=20x20cm

model UV2035-1 : 254nm, view size=20x35cm  
 model UV2035-2 : 306nm, view size=20x35cm  
 model UV2035-3 : 365nm, view size=20x35cm  
 model UV2035-B : 254 /365nm, view size=20x35cm

### Option accessories:

- OA-004-S UV transparent plate (25x25cm)
- OA-004-L UV transparent plate (25x45cm)
- OA-005-S UV transparent holder (7x9cm)
- OA-005-M UV transparent holder (10X12cm)
- OA-005-L UV transparent holder (12.5X15cm)
- OA-005-XL UV transparent holder (18X21cm)
- OA-006 UV transparent tube rack (5x4well)
- OA-007 Anti-UV goggle (eye protection)
- OA-008 Anti-UV face shield (face protection)
- OA-009 Gel cutting hood, black with window
- OA-010 UV transparent plate+holder(19x19cm)


UV block cover can fix at front or back side



Expand to image system compact size, save space



UV transparent plate, goggle



Gel holder



Gel cutting hood



萬柏科儀有限公司/666儀器網

www.lab666.com.tw

02-2451-4745

e-mail: sales@cae.url.tw